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CALIBRATION

Keeping pressed the **CALI** key for about 4 seconds, the "CAL" message will be displayed.

Press **SET** to display the weight calibrated last time. Press **↑** and **ENTER** to input the weight value you want to load (equal to at least 50% of the maximum quantity to be weighed).

Load the corresponding weight and press **ENTER** to confirm it.

Press **MODE** to exit from the calibration in any moment.

RS232 SERIAL PORT SETTINGS (OPTIONAL)

To enter serial port settings:

- turn off the instrument and then turn it back on
- while the display shows its test numbers form "999999" to "111111" simultaneously press the **TARE** and **MODE** keys
- the instrument beeps; release the keys
- at the end of the display test an acronym identifying the present serial operating mode appears

Now select the desired mode using the **↑** key:

- "CLOSE" to disable the serial port
- "-Stb-": stable weight transmission: the instrument start transmitting the string when the weight is stable and keeps on until the weight becomes unstable. (about 5 transmission per second)
- "-Etb-": press the **ENTER** key to send the string (only if the instrument is in weighing mode)
- "-Ser-": continuous weight transmission (about 5 transmission per second)

Confirm the selection by pressing **ENTER**: the selected baud rate will be displayed (1200, 2400, 4800, 9600). Press the **↑** key to change it, then press **ENTER** to confirm and exit the programming.

RS232 SERIAL PORT (DB9 CONNECTOR) PIN-OUT

Pin 2: TX

Pin 3: RX

Pin 5: Ground (GND)

The instrument can be connected to a PC by a pin-to-pin extension cable.

TRANSMISSION PROTOCOL

Serial line settings :

In order to enable the serial port, follow the instructions on the previous page.

The LACSN sends a formatted string containing to the value on the display at the speed of about 5 transmissions per second with these settings:

- 1 start bit
- 8 data bits
- 1 stop bit
- No parity

Transmitted string formatting:

The string contains 6 bytes, each byte has the following meaning:

Byte 0: String start character, hexadecimal FF

Byte 1: Each bit has the following meaning:

Bit 7	Bit 6	Bit 5	Bit 4	Bit 3	Bit 2	Bit 1	Bit 0
1= weight in Overflow	1= stable weight	1= negative weight	Operating mode 00 = weight 01 = count 10 = percentage		Decimal point position (number of decimals + 1) from 1 (no decimals) to 6 (5 decimals)		

Byte 2: less significant BCD part of the displayed value

Byte 3: intermediate significant BCD part of the displayed value

Byte 4: more significant BCD part of the displayed value

Byte 5: measurement unit (1= lb ; 0 = g or kg, depending on the instrument model).

In any case the transmitted value (weight, counting or percentage, depending on the instrument operating mode) is the one displayed and possible reset or tare are not taken into account.

Example:

Hexadecimal received string:

FF 42 55 08 00 00

- Byte 0: hex FF: string header
- Byte 1: hex 42:

Bit 7	Bit 6	Bit 5	Bit 4	Bit 3	Bit 2	Bit 1	Bit 0
0	1	0	0	0	0	1	0
No overflow	Stable weight	Positive weight	00: weight mode		Bin. 010 = dec. 2 → 1 decimal		

- Byte 2: hex 55: less significant part of the weight: 55
- Byte 3: hex 08: intermediate part of the weight : 08
- Byte 4: hex 00: more significant part of the weight: 0
- Byte 5: hex 00: measurement unit: g

Therefore the weight is: 000855; inserting the decimal point (1 decimal) → 85.5 g

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